


FLASH 31 Immobilization Kit



Complete immobilization system. By using together spinal board, spider straps and head immobilizer included in this kit, is possible to stabilize and immobilize the patient during rescue operations and ambulance transport.

Specific Characteristics

Lightweight and compact spinal board
Spider straps with strap closure
Head immobilizer with opening for ear inspection
Soft head immobilizer pillows to increase comfort
Wide strap surfaces make the head block placing easier
Adjustable frontal and chin straps with soft and high-adhesion paddings to optimize the immobilization
Smooth easy to clean surfaces

Kit components

ST02010B	Rock pin – Yellow spinal board with pin
ST02035A	Rock straps – Univ. spider straps for sp.board
SH00110A	Super blue – Universal head immobilizer

Accessories

ST02106A	Mark Up – system for permanent personalization
ST02105A	Fix board - 10 g vertical wall support for sp.board
ST02108A	Fix board - 10 g horiz. wall support for sp.board

Technical data

Spinal board dimensions	1840x450x50 ± 5mm
Spinal board weight	6,5 kg
Loading capacity	200 kg
Central belt length	1230 mm
Belt weight	460 g
Head block dimensions ¹	255x85x160 mm
Head immob. base dimensions ²	365x260x20 mm
Head immobilizer weight	1,46 ± 0,07 kg
Materials	PE, Nylon, PVC
¹ Size of one headblock	
² Buttonhole excluded	

Med.dev. CE marked according to dir. 93/42/CE

FLASH 31 ROCK PIN YELLOW +
ROCK STR. + SUPER BLUE

ST02048B



Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it