


## B-BAK PIN

### Spinal board with pins



B-bak Pin spine boards are devices designed to lift and immobilize patients with suspected spinal injuries.

#### Specific features

The number of handles, their placing and their shape offer an easy grip even in difficult conditions

The tapered shape allows to place the device inside Spencer basket stretchers

Smooth and rounded surfaces ensure an effective sanitation

Robust and light structure to ensure high resistance

Slides on the lower side to minimize sliding resistance during the application

The low thickness makes the application and storage easier

Equipped with pins to allow easy fixing of belts with snap hook

Allows the execution of preliminary X-ray diagnostic investigations

#### Technical data

Dimensions	1840 x 405 x 45 mm
Materials	PE
Number of handles	14
Number of pins	8
Load capacity	180 kg
Weight	6 kg

#### Accessories

ST02105A	Fix board – 10g vert. fastener for spine boards
ST02108A	Fix board – 10g horiz. fastener for spine boards
ST02035A	Rock Straps – Universal belt system
ST00494B	STX 494 – 3 belts with snap hook
ST02106A	Mark-Up – Permanent personalization system
SH00201A	Spencer Contour – universal head immobilizer

Class I MD compliant with UE Reg. 2017/745

B-BAK PIN YELLOW SPINE BOARD  
COMPLETE OF PINS



**ST02061B**

**CND Classification** V0804

**Registration number** 101249

Rev.0 (14/06/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222  
[www.spencer.it](http://www.spencer.it)