

Immobilization kit



Immobilization system complete of spinal board, head immobilizer and belts. This set allows to stabilize the patient during patient recovery procedures and ambulance transport.

Specific features

Extremely lightweight and compact spine board
Belts with strap closure
Head immobilizer with opening for ear inspection
Soft head immobilizer pillows to increase comfort
Wide strap surfaces make the head block placing easier
Adjustable frontal and chin straps with soft and high-adhesion paddings to optimize the immobilization
Smooth easy to clean surfaces

Technical data

Spinal board dimensions	1840x405x45 ± 5mm
Spinal board weight	6 ± 0,3 kg
Handles	14
Loading capacity	180 kg
Belts weight	460 g
Head block dimensions ¹	255x85x160 mm
Head immob. base dimensions ²	365x260x20 mm
Head immobilizer weight	1,5 ± 0,1 kg
Materials	PE, Nylon, PVC
¹ Size of one head block	
² Buttonhole excluded	

Accessories

ST02106A	Mark Up – Perm. Personalization system
ST02105A	Fix board - 10 g vertical wall support for sp.board
ST02108A	Fix board - 10 g horiz. wall support for sp.board

Med.dev. CE marked according to dir. 93/42/CE

B-BAK PIN W/SUPER BLUE AND
ROCK STRAPS

ST02065B



Rev.1 (27/04/21)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it