

## **B-BAK PIN MAX**

## Spinal board with pins



B-bak Pin Max spine boards are devices designed to lift and immobilize patients with suspected spinal injuries.

Specific features
The number of handles, their placing and their shape offer an
easy grip even in difficult conditions
The tapered shape allows to place the device inside Spencer
basket stretchers
Smooth and rounded surfaces ensure an effective sanitation
Robust and light structure to ensure high resistance
Slides on the lower side to minimize sliding resistance during the application
The low thickness makes the application and storage easier
Equipped with pins to allow easy fixing of belts with snap hook

Technical data		
Dimensions	1840 x 405 x 45 mm	
Materials	PE	
Number of handles	14	
Number of pins	8	
Load capacity	454 kg	
Weight	6 kg	

Accessories	
ST02105A	Fix board – 10g vert. fastener for spine boards
ST02108A	Fix board – 10g horiz. fastener for spine boards
ST02035A	Rock Straps – Universal belt system
ST00494B	STX 494 – 3 belts with snap hook
ST02106A	Mark-Up – Permanent personalization system
SH00201A	Spencer Contour – universal head immobilizer

Class I MD compliant with UE Reg. 2017/745		
B-BAK PIN MAX YELLOW SPINE		
BOARD	( (	
ST02067B		
CND Classification	V0804	