

Spine Pack

Immobilization kit



Immobilization kit including a practical carry bag. The spinal board, spider straps system and cervical collars included in this kit, allows to stabilize the patient during recovery procedures and ambulance transport.

Specific Characteristics

- Spine board characterized by high solidity
- Spider belts system to ensure a complete fixation of the patient
- Head immobilizer equipped with front and chin belts
- Head immobilizer with openings for auricular inspection
- Soft pillows to increase patient comfort
- Easy to clean thanks to smooth and filleted surfaces
- Selection of one-piece cervical collars Nec York
- Isothermic sheet suitable for hypothermia and burns

Technical data

| | |
|-----------------------------------------------|----------------------------|
| Spinal board dimensions | 1800x440x65 ± 10 mm |
| Spinal board weight | 7,4 ± 0,3 kg |
| Loading capacity | 192 kg |
| Belts weight | 460 g |
| Head immob.dimensions ¹ | 235x135x160 ± 5 mm |
| Head immob.weight | 1,5 kg |
| Materials | PE, PU, PVC |
| Included collar sizes | Paediatric, Short, Regular |
| ¹ Size referred to a single pillow | |

SPINE PACK - INTEGR. BAG
FOR SPINAL BOARD
COMPLETE
ST02103B



Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it