

## FIX BOARD

10G vertical wall support for spinal board



The image can differ from the actual product

Fastener designed for vertical fixation of Spencer spinal boards.

### Specific features

Extremely robust and lightweight  
Easy to use and install  
Quick locking and unlocking system

### Technical data

Upper part dimensions	570 x 25 x 65 ± 5 mm
Lower part dimensions	180 x 50 x 55 ± 5 mm
Materials	Steel, Nylon, Brass

### Reference standard

EN 1789:2007+A2:2014

Class I MD compliant with UE Reg. 2017/745

FIX BOARD 10 G WALL SUPPORT FOR  
SPINAL BOARD  
ST02105A



CND Classification	V9099
Registration number	99875

Rev.0 (14/06/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222  
[www.spencer.it](http://www.spencer.it)