

## FIX BOARD 10G

Wall support for spinal board



Fastener designed for horizontal fixation of Spencer spinal boards.

### Specific Characteristics

Extremely robust and lightweight  
Easy to use and install  
Quick locking and unlocking system

### Technical data

Dimensions <sup>(one piece)</sup>	545x25x65 ± 5 mm
Material	Steel, Nylon, brass
Weight <sup>(overall)</sup>	1,5 ± 0,2 kg

### Reference standard

EN 1789

Class I Med.Device according to Dir. 93/42/CEE

FIX BOARD 10G WALL SUPPORT FOR  
SPINAL BOARD

ST02108A

CND Classification



V9099

Registration number

806067

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)