

FIX BOARD

10G horizontal wall support for spinal board



Fastener designed for horizontal fixation of Spencer spinal boards.

Specific features

Extremely robust and lightweight
Easy to use and install
Quick locking and unlocking system

Technical data

Dimensions ^(one piece)	545 x 25 x 65 ± 5 mm
Materials	Steel, Nylon, Brass
Weight ^(overall)	1,5 ± 0,2 kg

Reference standard

EN 1789

Class I MD compliant with UE Reg. 2017/745

FIX BOARD 10G WALL SUPPORT FOR
SPINAL BOARD

ST02108A

CND Classification



V9099

Registration number

806067

Rev.0 (14/06/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it