

TANGO

Integrated pediatric and adult spine boards



Device for recovery, immobilization and transfer to the stretcher of adult and pediatric patients with suspected spinal injuries.

Specific features

The integration of adult and pediatric spinal boards, offers unmatched versatility

Pediatric spinal board with differentiated occipital areas for the best adaptation to the anatomical characteristics of pediatric patients

The pediatric spinal board can be easily extracted from the adult one, ensuring fast rescue operations and minimizing storage space

The head immobilizer can easily placed thanks to strap bands.

Technical data

Spinal board dimensions	1830x445x60 ± 5 mm
N° of handles	14
Pediatric sp.board weight	2 kg
Adult sp.board weight	8 kg
Weight of complete sp.board	10 kg
Pediatric sp.board dimensions	1195x320x45 ± 5 mm
Pediatric sp.board N° handles	10
Materials	PE, PVC, Nylon

Optionals

ST02102B	Spine pack – Carry bag
ST02106A	MARK UP – Permanent personalization system
ST02105A	FIX BOARD 10 G WALL SUPP.FOR SPINAL BOARD

Class I Med.Device according to Dir. 93/42/CEE

TANGO ADULT AND PAEDIATRIC SPINE BOARD



ST02140C

CND Classification	V0804
Registration number	101670

Rev.5 (02/11/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it