

FLASH 25

Immobilization system



Immobilization system intended for recovery and transfer on a main stretcher of adult and pediatric patients.

Specific Characteristics

The integrated adult and pediatric spine board, offers unmatched versatility.

Pediatric spine board with differentiated occipital areas to fit the anatomical characteristics of pediatric patients

Pediatric board equipped with specific belts with integrated meter to check the patient's height, allowing to identify which side of the board needs to be used

The base of Adult head immobilizer is padded to increase patient comfort

The head immobilizer fixing takes place by means of straps bands and is equipped with chin and forehead bands to improve head immobilization

Adult and pediatric head immobilizer have large side openings to allow ear inspection

The soft pillow padding increases patient comfort

Adult adjustable spider belts with strap closure provide greater level of immobilization and adaptability.

Technical data

Spinal board dimensions	1830x445x60 ± 5 mm
N° carry handles	14
Spinal board weight	9,2 kg
Pediatric spinal board	1195x320x45 ± 5 mm
N° carry handles pediatric board	10
Weight pediatric board	2,2 kg
Adult pillow dimensions ^(one piece)	235x135x160 ±10mm
Head immob. base dimensions	380x270x20 ±10mm
Adult Head immobilizer weight	1,5 ±0,1 kg
Pediatric pillow dimensions	190x70x120 ± 5mm
Pediatric pillows weight	560 ± 10g
Adult belts weight	750 ± 20 g
Pediatric belts weight	460 ± 10 g
Materials	PE, Nylon, TPU

Accessories

ST02102B	Spine pack – transport bag
ST02106A	MARK UP – Perm. Personaliz.system for sp. boards
ST02105A	FIX BOARD – 10G Vertical fastener for sp.boards

Med.dev. CE marked accordint to dir. 93/42/CEE

FLASH 25 - TANGO+ROCK STRAPS
AD/PED + CONTOUR+ PEDI GO

ST02147C



Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it