

FLASH 26

Immobilization system



Immobilization system intended for recovery and transfer on a main stretcher of adult and pediatric patients.

Specific Characteristics

The integrated adult and pediatric spin board, offers unmatched versatility.

Pediatric spine board with differentiated occipital areas to fit the anatomical characteristics of pediatric patients

Pediatric board equipped with specific belts with integrated meter to check the patient's height, allowing to identify which side of the board should be used

Head immobilizer and pediatric spinal board obtained by the adult one easily, allowing to be faster in rescue operations and to minimize storage space inside the ambulance

The head immobilizer fixing takes place by means of straps bands and is equipped with chin and forehead bands to improve head immobilization

Adult and pediatric head immobilizer have large side openings to allow ear inspection

The soft pillow padding increases patient comfort

Adjustable spider belts with strap closure provide greater level of immobilization and adaptability.

Accessories

ST02102B	Spine pack – transport bag
ST02106A	MARK UP – Perm. Personaliz.system for sp. boards
ST02105A	FIX BOARD – 10G Vertical fastener for sp.boards

Technical data

Spinal board dimensions	1830x445x60 ± 5 mm
N° carry handles	14
Spinal board weight	9,2 kg
Pediatric spinal board	1195x320x45 ± 5 mm
N° carry handles pediatric board	10
Weight pediatric board	2,2 kg
Head immobilizer dimensions	380x265x180 ± 10 mm
Head immobilizer weight	2070 ± 20 g
Adult belts weight	750 ± 20 g
Pediatric belts weight	460 ± 10 g
Materials	PE, Nylon, TPU

Med.dev. CE marked according to dir. 93/42/CEE

FLASH 26 - TANGO + TANGOFIX + REFLEX STRAPS + RSP

ST02148C



Registration number

43743

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it