

# FLASH 28

## Pediatric immobilization Kit



Pediatric immobilization system. By using together spinal board, spider straps and head immobilizer included in this kit, is possible to stabilize and immobilize the patient during rescue operations and ambulance transport.

### Specific features

- Pediatric spine board with differential occipital areas to adapt to the anatomical features of pediatric patients
- Spider straps with integrated meter to allow measurement of patient in order to decide which side of the spinal board should be used
- Head immobilizer fixed by straps and equipped with head and chin belts
- Head immobilizer with opening for ear inspection
- Head immobilizer with soft paddings to increase patient comfort
- Easy to clean thanks to smooth and rounded surfaces

### Accessories

- ST02144B – Wall support for Baby go

### Technical data

Spinal board dimensions	1195x320x45 mm
Spinal board weight	2,2 kg
Belts weight	350 g
Head block dimensions <sup>1</sup>	190x70x120 mm
Head immob. weight	650 g
Materials	PE, Nylon, PVC
Loading capacity	30 kg
<small><sup>1</sup>Size of one head block</small>	

FLASH 28 - BABY GO + PEDI  
GO + RSP STRAPS  
ST02190B



Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)