

Pedi Loc

Pediatric immobilization system



Pedi Loc spine boards are devices designed to lift and immobilize patients with suspected spinal injuries.

Specific features

- The number of bands and adjustment options, provide an effective immobilization
- The small size allows easy handling in tight spaces
- Removable cover for easy cleaning

Technical data

Dimensions	1220 x 250 x 30 mm
Number of carry handles	5
Head immob. dimensions	185 x 65 x 120 mm
Bands for	Shoulders, Thorax, Pelvis, Femur, Tibia
Materials	PVC, Wood, Polyurethane
Weight (with bag and accessories)	4,5 kg
Maximum load capacity	30 kg

Standard equipment

- Pedi Go – Pediatric head immobilizer
- Spine Pack – Carry bag

Class I MD compliant with UE Reg. 2017/745

PEDI LOC - PAEDIATRIC IMMOB. AND TRANSP. SYSTEM
ST02600C



Classification CND	V0804
Registration number	101701

Rev.0 (14/06/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it