

PEDI PACK

INTEGRATED SYSTEM FOR RESCUE
AND TRANSPORT OF TRAUMATIZED CHILDREN



Bag to store and transport the main devices commonly used in first aid procedures and recovery of traumatized children.

Specific characteristics

Possibility to insert a pediatric spinal board with a maximum size of 1250x300 mm

Compatible with Spencer spinal boards Baby-Go and Pedi Loc

The spinal board compartment, can be completely opened providing easy access to the content and has a zip closure.

Large pocket with rubberized mesh on the internal side of the front flap

3 separate compartments on the frontal side equipped with zip closure

The carrying handle is PVC covered for a better grip and to increase comfort during transport

The reinforced perimeter confers high resistance

Technical data

Main compartment dimensions	1250x300x150 mm
Frontal pocket dimension	290x210x110 mm
Overall volume	76 l
Materials	PVC, PP
Weight	1,61 ± 10g

PEDI PACK - PAEDIATRIC RESCUE / TRANSPORT SYSTEM
ST02700B

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it