


SHELL

Universal Basket stretcher



Spencer Shell is a basket stretcher intended for recovery of patients where other devices would be difficult to use. Typical uses are water rescue, recovery by ground hoisting, fast recoveries and transfers that require side impact protection.

Specific Characteristics

The single shell body makes it extremely safe and easy to handle

The perimetral frame provides rigidity and allows to have a firm grip on the handles

Equipped with rings for hooks fixing

Perimetral rope to fix belts or accessories

The shape of the inferior part reduces friction during ground dragging procedures

Mattress to increase patient comfort

Equipped with straps and rings

Accessories

ST04525B	STX 537 – One piece universal buoyant
ST04518A	STX 518 – Two piece universal buoyant
ST04524B	STX 538 – Three piece universal buoyant
ST04519C	STX 519 – Adjustable lifting bridle
ST04522B	STX 540 – Fixed lifting bridles
ST04523B	STX 542 – Pulling bridles
ST04600B	STX 580 – Trolley for basket stretchers
ST04601B	STX 581 – Trolley for basket stretchers w/big wheels
ST04040A	STX 40 – Carry bag for Shell
IF01059B	Track 5 – IV holder for basket stretchers

Standard Equipment

ST04011A	Footrest
ST00598A	STX 598 – RESTRAINT BELTS

Technical data

Length	2150 ± 10 mm
Width	650 ± 10 mm
Height	200 ± 10 mm
N° of handles	12
Materials	PE, PP, Al
Weight ^(without standard equipment)	12,5 ± 0,5 kg
Loading capacity	280 kg

Class I Med.Device according to Dir. 93/42/CEE

SPENCER SHELL BASKET
STRETCHER MILITARY



ST04001A

CND classification V0899

Rev.0 (10/10/2019)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it