

FOOTREST

For basket stretchers



Footrest intended to limit the longitudinal movements of the patient during the use of Spencer basket stretchers.

Specific Characteristics

Lightweight and durable body with large bearing surface

Easy to clean

Flexibility of use guaranteed by the additional hook

Technical data

Shell dimensions	350x65x160 ± 10 mm
Distance between hooks	Min 560 mm Max 660 mm
Material	PE, Nylon, Steel
Weight	420 ± 20g

Class I Med.Device according to Dir. 93/42/CEE

FOOTREST FOR SHELL AND
TWIN SHELL

ST04011A

CND Classification



V0899

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it