

STX 40

Spencer Shell carry bag



Bag for transport and storage of Spencer Shell basket stretchers.

Specific features

PVC fabric easy to clean
 Handles coated with soft plastic material for easy transport
 Longitudinal belts with ring terminal allow wall hanging
 Upper strap closure

Technical data

Dimensions	2150 x 610 x 200 ± 10 mm
Number of handles	6
Material	PVC
Weight	2,35 ± 0,05 kg

STX 40 - BAG FOR SHELL IN
 YELLOW PVC
 ST04040A

Rev.0 (27/11/2020)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it