

SPENCER SHELL

Universal basket stretcher



Spencer Shell basket stretchers are intended for the recovery and transport of patients. They can be used for rescue operations in those cases where it is necessary to protect patients from side impacts, where it is not possible to intervene with normal transport equipment. The devices can be hoisted with lifting systems applied to fixed anchors, in all cases keeping the stretcher in a horizontal position with respect to the ground, and only if equipped with a Spencer harnessing system. It is not foreseen that the patient be able to intervene on the devices.

Specific features

The single shell body makes it extremely safe and easy to handle

The perimetral frame provides rigidity and allows to have a firm grip on the handles

Equipped with rings for hooks fixing

Perimetral rope to fix belts or accessories

The shape of the inferior part reduces friction during ground dragging procedures

PVC mattress to increase patient comfort and to allow easy cleaning

Accessories

ST04525B	STX 537 – One piece universal buoyant
ST04518A	STX 518 – Two piece universal buoyant
ST04519C	STX 519 – Adjustable lifting bridles
ST04522B	STX 540 – Fixed lifting bridles
ST04040A	STX 40 – Carry bag for Spencer Shell

Standard equipment

STX 592 – Restraint belts 2 pcs with metal buckle (three belts are supplied)

Footrest

Technical data

Lenght	2150 ± 10 mm
Width	650 ± 10 mm
Height	200 ± 10 mm
Number of handles	12
Materials	PE, PP, Al
Weight ^(without accessories)	12,5 ± 0,5 kg
Load capacity	280 kg

Class I MD compliant with UE Reg. 2017/745

MOD.S003 SHELL ORANGE,3
BELTS,FEETREST,YELLOW MAT.

ST04052

CND Classification



V0899

Rev.0 (27/11/2020)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it