

DAKAR

Basket stretcher with wheels



Spencer Dakar basket stretchers are intended for the recovery and transport of patients. They can be used for rescue operations in those cases where it is necessary to protect patients from side impacts, where it is not possible to intervene with normal transport equipment. The devices can be hoisted with lifting systems applied to fixed anchors, in all cases keeping the stretcher in a horizontal position with respect to the ground, and only if equipped with a Spencer harnessing system. It is not foreseen that the patient be able to intervene on the devices.

Specific features The single shell body makes it extremely safe and easy to handle Aluminum frame with telescopic handles increases safety during lifting operations and make them easier Equipped with eyelets for the attachment of snap hooks Perimetral rope to fix belts or accessories Wheels of large diameter on the front side of the stretcher to allow short ground displacements Mattress to increase patient comfort Footrest to improve patient's longitudinal containment

Accessories		
ST04519C	STX 519 – Adjustable lifting bridles	
ST04522B	STX 540 – Fixed lifting bridles	

	Standard equipment
	STX 592 – Restraint belts 2 pcs with metal buckle (three
	belts are supplied)
	Footrest

Technical data	
Length	2240 mm
Length with handles opened	2760 mm
Width	640 mm
Height	190 mm
Wheels	2 x ø 200 mm
Number of handles	12
Materials	PE, PP, Al
Weight (without accessories)	16 ± 0,5 kg
Load capacity	356 kg

Class I MD compliant with UE Reg. 201	7/745
DAKAR BASKET STRETCHER WITH	
WHEELS ORANGE	((
ST04070B	
CND Classification	V0899
Registration number	171997