

TOTAL

Recovery stretcher



Total recovery stretchers are devices to be used when the morphological characteristics of the intervention area require particularly compact equipment and adherence to the patient.

Total can also be lifted with the aid of hoisting systems fixed to the ground, using the dedicated harness system.

Specific features

- Belts for towing and lifting
- A system of belts for the assembling of the stretcher and subsequently folding before replacing it in the bag
- A system of additional handles for the transport of heavy patients by four or eight operators
- Aluminum hook to allow for lifting
- Transport bag in nylon Ripstop
- Hoisting harness in nylon

Technical data

Stretcher dimensions	2430 x 920 x 3 ± 10 mm
Bag dimensions	∅ 380 x 950 ± 10 mm
Stretcher weight	6 ± 0,5 kg
Weight complete system	10 ± 1 kg
Loading capacity	200 kg
Stretcher material	Polyethylene

Class I MD compliant with UE Reg. 2017/745

TOTAL YELLOW - RECOVERY
STRETCHER
ST04090A



CND Classification	V9099
Registration number	194201

Rev.0 (24/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it