

TOTAL BAG

Backpack carry bag for Total



Robust bag for transport and storage of Total recovery stretcher.

Specific characteristics

- Allows to reduce storage space
- Easy to carry thanks to the padded shoulder straps
- Pocket for storage of accessories
- Reflective label increases visibility
- Easy to clean

Technical data

Dimensions	ø380x950 mm
Materials	PVC, Nylon
Weight	1,5 kg

TOTAL BAG – BACKPACK CARRY BAG FOR
TOTAL

ST04093C

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it