

TOTAL KON

Recovery stretcher



Total Kon recovery stretchers are devices to be used when the morphological characteristics of the intervention area require particularly compact equipment and adherence to the patient.

Specific features

- Allows for reduced storage space
- Handles facilitate transportation
- Equipped with patient handles
- Easy to sanitise

Technical details

Material	Polyethylene
Length	1200 mm
Width	918 mm
Thickness	3 mm
Loading weight	200 kg
Weight	2 kg
Transport handles	Polyester
Rope for prolonged towing	Nylon
Harness bands for the patient	Polyester
Colour	Yellow

Class I MD compliant with UE Reg. 2017/745

TOTAL KON - HALF BODY RESCUE
STRETCHER
ST04095A



Rev.0 (24/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it