

STX 10

Twin Shell carry bag



Bag for transport and storage of Spencer Twin Shell basket stretchers.

Specific features

PVC fabric easy to clean

Handle coated with soft plastic material for easy transport

Longitudinal belts with ring terminal allow wall hanging

Upper strap closure

Technical data

Dimensions	1200 x 630 x 280 ± 10 mm
Material	PVC
Weight	1,64 ± 0,04 kg

STX 10 - YELLOW PVC BAG FOR
TWIN SHELL
ST04110A

Rev.0 (27/11/2020)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it