

## BOSTON TEC

### Stainless steel basket stretcher



Boston TEC is a basket stretcher intended for the recovery of obese patients in cases where other devices would be difficult to use. Typical uses are recoveries with the help of winches on the ground, recoveries and transfers that require rapid intervention and protection from side impacts.

#### Specific features

The stainless steel structure makes it particularly robust and resistant

The shape of the frame provides a firm grip to the operators together with a wide choice of grip points

The polyethylene board increases thermal insulation and allows easy sanitation

Trendelenburg antishock position

Adjustable backrest

#### Standard equipment

ST70002A | 3 belts two pieces with metal buckle

#### Technical data

Length	2110 ± 10 mm
Width	650 ± 10 mm
Height	185 ± 10 mm
Materials	Steel, PE, Nylon
Weight <small>(without standard equipment)</small>	23 ± 0,5 kg
Loading capacity	360 kg

#### Optional accessories

ST04519C	STX 519 – Adjustable harness
ST04522B	STX 540 – Standard harness
ST04523B	STX 542 – Towing belts

Class I Med. Dev. according to Dir. 93/42/CEE

BOSTON TEC STAINLESS STEEL  
BASKET STRETCHER YELLOW  
**ST04300B**



CND Classification	V0899
Registration number	172055

Rev.0 (10/10/2019)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)