BOSTON PRO Basket stretcher with spinal board



Spencer Boston Pro basket stretchers are intended for the recovery and transport of patients. They can be used for rescue operations in those cases where it is necessary to protect patients from side impacts, where it is not possible to intervene with normal transport equipment. The devices can be hoisted with lifting systems applied to fixed anchors, in all cases keeping the stretcher in a horizontal position with respect to the ground, and only if equipped with a Spencer harnessing system. It is not foreseen that the patient be able to intervene on the devices.

Specific features

- The aluminum structure makes it particularly robust, resistant and reduces the weight
- The shape of the frame provides a firm grip to the
- operators together with a wide choice of grip points
- The removable spinal board makes it very versatile in
- rescue operations

Accessories	
ST04519C	STX 519 – Adjustable lifting bridles
ST04522B	STX 540 – Fixed lifting bridles

Standard equipment

STX 702 Three reflective belts, two pieces, with metal buckle

Technical data	
Length	2110 mm
Width	650 mm
Height	250 mm
Materials	Al, PE, Nylon
Weight (without accessories)	17 kg
Load capacity	360 kg

Class I MD compliant with UE Reg. 2017/745	
BOSTON PRO ALUMINIUM BASKET	
STRETCHER WITH SPINAL BOARD	
ST04302B	
CND Classification	V0899

Rev.0 (27/11/2020) UNCHECKED COPY – further revisions will be available on <u>http://support.spencer.it</u> Spencer Italia S.r.I. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it