

BOSTON PRO

Basket stretcher with spinal board



Boston PRO is a basket stretcher intended for the recovery of obese patients in cases where other devices would be difficult to use. Typical uses are recoveries with the help of winches on the ground, recoveries and transfers that require rapid intervention and protection from side impacts.

Specific features

The stainless steel structure makes it particularly robust and resistant

The shape of the frame provides a firm grip to the operators together with a wide choice of grip points

The removable spinal board makes it very versatile in rescue operations

Standard equipment

ST70002A | 3 belts two pieces with metal buckle

Technical data

Length	2110 ± 10 mm
Width	650 ± 10 mm
Height	250 ± 10 mm
Materials	Steel, PE, Nylon
Weight <small>(without standard equipment)</small>	26 ± 0,5 kg
Loading capacity	360 kg

Optional accessories

ST04519C	STX 519 – Adjustable harness
ST04522B	STX 540 – Standard harness
ST04523B	STX 542 – Towing belts

Class I Med. Dev. according to Dir. 93/42/CEE

BOSTON PRO BASKET STRETCHER
WITH SPINAL BOARD
ST04303B



CND Classification

V0899

Rev.1 (10/10/2019)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it