


STX 518Universal floating devices two pieces for basket stretchers or spine boards





The floating devices are accessories for basket stretchers or spine boards, to be used to increase the degree of buoyancy in water of the device with which they are used. The device is not suitable for use with the Dakar, Dakota and Dakota Light basket stretchers because they cannot be used in water.

Specific features
Fastened by means of 3 quick release buckles
PVC removable coating for easier cleaning and drying

Technical data	
Floatation device dimensions	\emptyset 150 x 1020 ± 10 mm (per piece)
Total buoyancy thrust With floatation device fully immersed	300 N
Materials	PE, PVC, Nylon
Weight	$2,2 \pm 0,1 \text{ kg}$

Class I MD compliant with UE Reg. 20	17/745
STX 518 2PIECES-UNI.FLOAT.SYS.FOR	
BASK. STRETCH.	(6
ST04518B	
CND Classification	V0899
Registration number	188130