STX 519 Adjustable lifting bridles system



These harnesses are accessories for basket stretchers, to be used solely for lifting from a floor-mounted position (STX 519 and STX540) or for dragging (STX 519). It is not foreseen that the patient be able to intervene on the devices.

Specific features	Technical data	Technical data	
4 hook with automatic locking secure Easy to adjust system	Belt length	Minimum 1250 mm Maximum 1680 mm	
	Belt width	50 mm / 20 mm	
	Weight	1.8 ± 0.1 kg	
	Material	Nylon, Steel , Aluminum	
	Fastening type	Spring catch	
	Number of fastening bands	4	
	Compatibility	All basket stretchers manufactured by Spencer	
	Maximum applicable static safety load	450 kg (1000 lbs)	
		Class I MD compliant with UE Reg. 2017/745 STX 519 ADJUST. LIFTING BRID.	
	- LIFTING BRIDLE SYS		

CND ClassificationV0899Registration number188142

Rev.0 (27/11/2020) UNCHECKED COPY – further revisions will be available on <u>http://support.spencer.it</u> Spencer Italia S.r.I. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it