


STX 537

Universal floating device for basket stretchers or spine boards



The floating devices are accessories for basket stretchers or spine boards, to be used to increase the degree of buoyancy in water of the device with which they are used. The device is not suitable for use with the Dakar, Dakota and Dakota Light basket stretchers because they cannot be used in water.

Specific features

Float with fastening by means of quick release buckles
Removable PVC cover to make cleaning and drying easier

Technical data

Floatation device dimensions	Ø 150 x 3330 ± 10 mm
Total buoyancy thrust With floatation device fully immersed	550 N
Materials	PE, PVC, Nylon
Weight	3,1 ± 0,1 kg

Class I MD compliant with UE Reg. 2017/745

STX 537 1PIECE-UNI. FLOAT. SYST. FOR
BASK. STRETCH
ST04525B



CND Classification	V0899
Registration number	187804

Rev.0 (27/11/2020)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it