

SPENCER SX-R

Pickup stretcher



Spencer SX is a pickup stretcher with two blades for lifting and loading of patients on immobilising or transport devices such as spinal boards, vacuum mattresses or transport stretchers.

Specific characteristics

- Extremely light and resistant structure
- The possibility to adjust the length confers high versatility
- Foldable to save space when stored
- Smooth surfaces easy to clean
- The numerous handles make the transport easier
- Belts with transport bag

Standard accessories

- | | |
|----------|---------------------------------------|
| ST00593A | K belt 3—three one piece belts STX597 |
|----------|---------------------------------------|

Optional accessories

- | | |
|----------|------------------------------------|
| ST05009A | SX Bag – Orange transport bag |
| ST05018C | SXV – 20G vertical fixing system |
| ST05019C | SXO – 20G horizontal fixing system |

Compliant with

UNI EN 1865-1

Technical data ¹

Shortened length	1660 mm
Intermediate length 1	1780 mm
Intermediate length 2	1900 mm
Intermediate length 3	2050 mm
Folded length	1190 mm
Thickness folded	70 mm
Thickness opened	60 mm
Width	440 mm
Material	Aluminum
Loading capacity	170 kg
Weight	7,6 ± 0,4 kg

¹ Length measures are subjected to tolerances of ± 10 mm

Class I Med.Device according to Dir. 93/42/CEE

SPENCER SXR PICKUP STR.2

BLADES W/SET 3 BELTS

ST05002A

CND Classification V9099

Registration number 91818



Rev.3 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it