

PICKUP STRETCHER



Pickup stretcher with two blades for lifting and loading of patients on immobilizing or transport devices such as spinal boards, vacuum mattresses or transport stretchers.

Specific features

- Extremely light and resistant structure
- hinged and joined arrangement, which is easy to set the stretcher into two parts, is set in the middle at both ends of the stretcher
- Adjustable length according to patients' body.
- Smooth surfaces easy to clean
- The numerous handles make the transport easier

Technical data ¹

| | |
|--------------------|-------------|
| Dimensions 1 (max) | 204x42x6 cm |
| Dimensions 2 | 192x43x6 cm |
| Dimensions 3 | 180x42x6 cm |
| Dimensions 4 (min) | 167x42x6 cm |
| Dimensions folded | 120x42x8 cm |
| Material | Aluminum |
| Loading capacity | 159 kg |
| Weight | 7 ± 0,4 kg |

¹ Length measures are subjected to tolerances of ± 10 mm

Class I Med.Device according to Dir. 93/42/CEE

PICKUP STRETCHER

ST05005D

Rev.0 (10/04/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
spencer.it