

SXV

Vertical fastener for Pick Up stretcher



Fastener studied to fix Spencer SX Pick Up stretchers in vertical position.

Specific Characteristics

Extremely robust and lightweight
Easy to use and install
Quick locking and unlocking system

Technical data

Dimensions ^(long piece)	440x25x95 ± 5 mm
Dimensions ^(short piece)	180x60x60 ± 5 mm
Material	Al, Nylon, brass
Weight ^(overall)	520 ± 20 g

Reference standard

NA

Class I Med.Device according to Dir. 93/42/CEE

SXV VERTICAL BRACKETS FOR
SPENCER SX

ST05017A

CND Classification



V9099

Registration number

91873

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it