

## SXO 20G

Horizontal fastener for Pick Up stretcher



Fastener studied for a safe fixing of Spencer SX Pick Up stretchers.

### Specific Characteristics

Extremely robust and lightweight  
Easy to use and install  
Designed to withstand accelerations up to 20 g

### Technical data

Dimensions	510x25x60 ± 5 mm
Material	Steel, Nylon
Weight	1,5 ± 0,2 kg

### Compliant with

EN 1789

Class I Med.Device according to Dir. 93/42/CEE

SXO HORIZONTAL WALL BRACK. 20G  
FOR PICK UP STRETCH

ST05019C

CND Classification V9099

Registration number 91927



Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)