

SX BAG

Pickup stretcher bag



Bag for transport and storage of Spencer Pickup stretchers.

Specific Characteristics

- PVC fabric easy to clean
- Handle coated with soft plastic material for easy transport
- Longitudinal belts with ring terminal allow wall hanging
- Upper strap closure

Technical data

Dimensions	1220x390x95 ± 10 mm
N° of handles	1
Material	PVC
Weight	960 ± 20g

TRANSPORT BAG FOR PICK UP STRETCHER
ST05055A

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it