

Jack S

Hydraulic multilevel support frame for stretchers



Jack S is intended for the handling and support of stretchers. It reduces the effort linked to the lifting of loads.

Specific Characteristics

Hydraulic lifting system
 Easy height adjustment thanks to the pedals located in easily accessible area
 2 fixed wheels and 2 pivoting with large diameter to make movement easier
 Polished stainless steel frame to facilitate cleaning operations
 2 brakes
 High load capacity

Technical data

Length	1900 mm
Width	630 mm
Weight	28 kg
Height	From 465 to 945 mm
Front wheels	∅ 200 mm
Rear wheels	∅ 260 mm
Load capacity	270 kg

Class I Med.Device according to Dir. 93/42/CEE

JACK S MULTILEVEL SUPPORT
 FRAME FOR STRETCHERS
ST10260B



CND Classification	V9099
Registration number	421142

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it