

Spencer 420

Transport chair with two wheels



Spencer 420 is a transport chair for the evacuation and transport of patient from the rescue place to the ambulance or to another device intended for transport on the vehicle.

Specific Characteristics

Easy to clean PVC sheets
 2 rear foldable handles to reduce encumbrance and to facilitate the lifting in every situation
 Telescopic front handles to facilitate the handling
 Foldable to save space
 Strong and light frame

Standard equipment

ST00423A Pair of belts for chair 420/425 orange

Accessories

ST30451A MOUNTING BRACKETS FOR CHAIR SPENCER 450-55/420-25
 ST36452E UPPER FIXATION SYSTEM FOR TRANSPORT CHAIRS

Technical data

Length opened w/closed handles	660 ± 10 mm
Length opened w/opened handles	1150 ± 10 mm
Width	510 ± 10 mm
Height opened	900 ± 10 mm
Thickness folded	210 ± 10 mm
Rear wheels	ø125 mm
Materials	Steel/Al/PU/PVC
Weight	10,5 kg ± 0,3
Loading capacity	180 kg

Class I Med.Device according to Dir. 93/42/CEE

STAIR CHAIR SILVER SPENCER 420S

ST10420A

CND Classification V9099

Registration number 93511



Rev.3 (02/04/2020)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it