

Spencer 425

Transport chair with four wheels



Spencer 425 is a transport chair for the evacuation and transport of patient from the rescue place to the ambulance or to another device intended for transport on the vehicle.

Caratteristiche specifiche

Easy to clean PVC sheets
2 rear foldable handles to reduce encumbrance and to facilitate the lifting in every situation
Telescopic front handles to facilitate the handling
Foldable to save space
Strong and light frame
Foldable footrest to increase stability and comfort of the patient
Front wheels with brake for rotation and axis

Standard equipment

ST00427A	Pair of belts for chair 420/425
----------	---------------------------------

Accessories

ST30451B	MOUNTING BRACKETS FOR CHAIR SPENCER 450-55/420-25
ST36452E	UPPER FIXATION SYSTEM FOR TRANSPORT CHAIRS
ST00428B	Stabilizer stuffed seat/backrest weld series

Technical data

Length opened w/closed handles	700 ± 10 mm
Length opened w/opened handles	1150 ± 10 mm
Width	490 ± 10 mm
Height opened	900 ± 10 mm
Thickness folded	210 ± 10 mm
Rear wheels	ø125 mm
Front wheels	ø100 mm
Materials	Steel/Al/PU/PVC
Weight	11 kg ± 0,3
Loading capacity	150 kg

Class I Med.Device according to Dir. 93/42/CEE

SPENCER 425 SILVER FOLDING
CHAIR 4 WHEELS ORANGE C
ST10425B



CND Classification	V9099
Registration number	93530

Rev.3 (28/01/2020)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it