

## **Spencer 4BELL Silver**

## Transport chair with four wheels





Spencer 4BELL Silver is a transport chair for the evacuation and transport of patient from the rescue place to the ambulance or to another device intended for transport on the vehicle.

Specific features
Easy to clean PVC sheets
2 rear foldable handles to save space and to facilitate the lifting in every situation
The telescopic handles locking system can be unlocked by a simple directional pressure
Foldable to save space
Strong and light frame
Foldable footrest to increase stability and comfort of the patient
Self-positioning front wheels with brake for rotation and axis
Soft backrest to increase patient comfort

Technical data	
Length w/ handles closed	750 ± 5 mm
Length w/ handles opened	1170 ± 5 mm
Width	485 ± 5 mm
Height opened	905 ± 5 mm
Dimensions folded	310x485x915 ± 5 mm
Rear wheels	ø150 mm
Front wheels	ø100 mm
Materials	Al/Steel/PU/PVC
Weight	9,4 kg ± 0,2
Loading capacity	180 kg

## Standard equipment Integrated belts

Accessories
ST10432B Fastener 4BELL MAX 10g certified

Reference standards

UNI EN 1865-1

UNI EN 1789 if used with the dedicated fastener

Class I Med.Device according to Dir. 93/42/CEE
TRANSPORT CHAIR 4BELL SILVER
ST10431B
CND Classification
V9099