

4BELL MAX

10G fastener for chair 4BELL



Fastener intended for installation inside an ambulance, to be used with 4BELL chair. The device can be fixed both on wall or floor, according to the vehicle structure and equipment. The device is a standard equipment of the chair 4BELL.

Specific Characteristics

The polyethylene body, dampens vibrations and makes cleaning procedures easier

Its special shape allows to fix safely also 4BELL chairs equipped with optional headrest and avoids incorrect placing of the chair

The chair is easily fixed simply by turning the red knob

Compliant with

EN 1789

Technical data

Dimensions	505x450x220 ± 5mm
Materials	PE, Nylon, Al
Weight	1,8 ± 0,2 kg

Class I Med.Device according to Dir. 93/42/CEE

FIXATION SYSTEM 4BELL MAX 10G

TESTED

ST10432B

CND Classification



V9099

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it