

## Spencer 4BELL COMFORT

### Transport chair with four wheels



Spencer 4BELL COMFORT is a transport chair for the evacuation and transport of patient from the rescue place to the ambulance or to another device intended for transport on the vehicle.

#### Specific features

- Easy to clean PVC sheets
- 2 rear foldable handles to save space and to facilitate the lifting in every situation
- Front telescopic handles make lifting easier
- The telescopic handles locking system can be unlocked by a simple directional pressure
- Foldable to save space
- Strong and light frame
- Foldable footrest to increase stability and comfort of the patient
- Self-positioning front wheels
- Soft backrest to increase patient comfort
- Independent brakes on rear wheels are easy to operate
- Foldable headrest to increase patient comfort

#### Technical data

Length w/ handles closed	750 ± 5 mm
Length w/ handles opened	1170 ± 5 mm
Width	510 ± 5 mm
Height with headrest closed	910 ± 5 mm
Dimensions folded	930x510x310 ± 5 mm
Seat height	480 mm
Seat dimensions	390x390 mm
Backrest height from the seat	460 mm
Backrest width	440 mm
Rear wheels	ø150 mm
Front wheels	ø100 mm
Materials	Al/Steel/PU/PVC
Weight <sup>(with standard equipment)</sup>	12 ± 0,2 kg
Loading capacity	180 kg

#### Standard equipment

- Integrated belts
- Headrest

Class I Med.Device according to Dir. 93/42/CEE

TRANSPORT CHAIR 4BELL COMFORT

ST10436

CND Classification



V9099

#### Accessories

- Fastener 4BELL MAX 10g certified

#### Reference standards

- UNI EN 1789 if used with the dedicated fastener

Rev.1 (18/07/2017)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)