

4BELL STAIR MAX

10G fastener for chair 4BELL STAIR



Fastener intended for installation inside an ambulance, to be used with 4BELL STAIR chair. The device should be fixed on wall, according to the vehicle structure and equipment.

Specific features

- The aluminum structure confers high resistance and a reduced weight
- Easy to install
- The chair is easily fixed simply by turning the red knob

Technical data

Dimensions	495x340x225 ± 5mm
Materials	Nylon, Al
Weight	1,33 ± 0,1 kg

Compliant with

EN 1789

Class I Med.Device according to Dir. 93/42/CEE

FASTENER 4BELL STAIR MAX 10G

ST10441B



CND Classification V9099

Registration number 1530324

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it