

4BELL+ MAX

10G fastener for chair 4BELL +





Fastener intended for installation inside an ambulance, to be used with 4BELL+ chair. The device can be fixed both on wall or floor, according to the vehicle structure and equipment. The device is a standard equipment of the chair 4BELL+.

Specific Characteristics
The single piece body allow easy installation and cleaning
The chair is easily fixed simply by turning the red knob

Fechnical data			
Dimensions	490x143x780 ± 5mm		
Materials	Al, PP, steel		

Compliant with		
	EN 1789	

Class I Med.Device according to Reg.	. 2017/745/UE
4BELL+MAX FASTENING SYSTEM	
ST10453B	CE
CND Classification	V9099