

Spencer 480

Comfort chair with four wheels



Spencer 480 is a transport chair for the evacuation and transport of patient from the rescue place to the ambulance or to another device intended for transport on the vehicle.

Specific features

- PVC sheet easy to clean and padded to maximize patient comfort
- 2 rear foldable handles to reduce overall dimensions and facilitate the lifting in every situation
- Telescopic front handles make lifting easier
- Foldable to save storage space
- Robust and lightweight structure
- Footrest to improve stability and patient comfort
- Swiveling front wheels equipped with brake for rotation and axis
- The large diameter of the rear wheels allow to absorb roughness of the terrain minimizing stress on the patient's back
- Padded armrests to improve containment and comfort

Technical data

Length w/ handles closed	785 ± 10 mm
Length w/ handles opened	1245 ± 10 mm
Width	565 ± 10 mm
Height	1175 ± 10 mm
Thickness folded	350 ± 10 mm
Rear wheels	ø200 mm
Front wheels	ø125 mm
Materials	Al/Steel/Nylon/PU/PVC
Weight	14 kg ± 0,3
Loading capacity	150 kg

Reference standard

NA

Standard equipment

ST00427A BLACK PAIR OF BELTS

Spare parts

ST10481C PVC SHEETS FOR CHAIR 480

Class I Med.Device according to Dir. 93/42/CEE

SPENCER 480 FOLDABLE CHAIR W/
ARMRESTS,WHEELS 200,
ST10480B



CND Classification V9099

Registration number 93312

Rev.5 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it