

Spencer 4BELL Silver + Transport chair with four wheels



Spencer 4BELL Silver is a transport chair for the evacuation and transport of patient from the rescue place to the ambulance or to another device intended for transport on the vehicle.

Specific features

- Easy to clean PVC sheets
- 2 rear foldable handles to save space and to facilitate the lifting in every situation. The locking system of the handles in opened position, increase safety during transport
- The telescopic handles locking system can be unlocked by a simple directional pressure
- Foldable to save space
- Strong and light frame
- Foldable footrest to increase stability and comfort of the patient
- Pivoting front wheels, mounted on a reinforced frame that makes the device even more maneuverable
- Soft backrest to increase patient comfort
- Easy to identify and easy to operate rear brakes

Standard equipment

- Integrated belts

Accessories

- ST10453 – 4BELL+ MAX 10G FASTENER
- ST10452 – SEAT XL
- ST10450 - Armrests
- ST10451 - Rear auxiliary handles

Technical data

Length w/ handles closed	760 ± 5 mm
Height opened	905 ± 5 mm
Length w/ handles opened	1120 ± 5 mm
Width	490 ± 5 mm
Dimensions folded	290x490x905 ± 5 mm
Rear wheels	Ø150 mm
Front wheels	Ø100 mm
Materials	Al/Steel/PU/PVC
Weight	9,4 kg ± 0,2
Loading capacity	180 kg

Reference standards

- UNI EN 1865-1
- UNI EN 1789 if used with the dedicated fastener

Class I Med.Device according to Regulation 2017/745/UE

4BELL SILVER+

ST10500

CND Classification



V9099

Rev.0 (14/03/2022)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it