

Spencer 4BELL Avio +

Transport chair with four wheels



Spencer 4BELL Avio+ is a transport chair for the evacuation and transport of patient from the rescue place to the ambulance or to another device intended for transport on the vehicle.

Specific features

- Easy to clean PVC sheets
- 2 rear foldable handles to save space and to facilitate the lifting in every situation. The locking system of the handles in opened position, increase safety during transport
- The telescopic handles locking system can be unlocked by a simple directional pressure
- Foldable to save space
- Strong and light frame
- Foldable footrest to increase stability and comfort of the patient
- Pivoting front wheels, mounted on a reinforced frame that makes the device even more maneuverable
- Soft backrest to increase patient comfort
- Easy to identify and easy to operate rear brakes

Standard equipment

Integrated belts

Accessories

ST10450	Armrests
ST10451	Rear auxiliary handles

Technical data

Length w/ handles closed	750 ± 5 mm
Height opened	905 ± 5 mm
Length w/ handles opened	1170 ± 5 mm
Width	460 ± 5 mm
Dimensions folded	310x930x460 ± 5 mm
Rear wheels	ø150 mm
Front wheels	ø100 mm
Materials	Al/Steel/PU/PVC
Weight	9,3 kg ± 0,2
Loading capacity	180 kg

Class I Med.Device according to Regulation 2017/745/UE

4BELL AVIO+

ST10502

CND Classification



V9099

Rev.0 (18/07/2022)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it