

## ELLA SELF

### Self-loading stretcher



Ella series field stretchers Ella Self are devices for lifting and carrying people for use in field rescue situations. The foldable, lockable legs allow for transferring either by means of wheels on the legs or by lifting when the legs are locked in the closed position. Parking is permitted with the stretcher alone or with a trolley with closed legs.

#### Specific features

Polyethylene board easy to clean
Color coded commands
Adjustable backrest in 5 positions in addition to the horizontal
Board with Trendelenburg antishock and Fowler position
2 fixed wheels and 2 pivoting with large diameter to make movement easier
Brakes on rear wheels
Foldable side bars <i>Safe bar</i>
System to lock the legs in closed position to lift and carry the device without using the wheels

#### Accessories

EN90003C	END-T Table for stretchers 10 g certified
ST00491A	STX 90 Telescopic headrest for stretchers
IF01044B	Track 4 telescopic IV holder for Ella Self
CB09025C	Tanker oxygen tank holder for stretcher side bar
CB09028C	Tanker support for Safe Bar
ST70000A	QMX 777 Black PVC anatomic mattress
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Thorax Belt w/integrated re-winding system
ST42100A	FP Posterior fastener for self-loading stretchers
ST42200A	FA Anterior fastener for self-loading stretchers

#### Standard equipment

STX 592 – Two piece restraints with metal buckle

#### Technical data

Length with handles closed	1980 ±10 mm
Width	570 ± 10 mm
Wheels	Ø 200 mm
Dimensions sidebars	680 x 200 mm
Materials	Steel/Nylon/Al/PE
Weight	40 ± 1 kg
Capacity	170 kg
Loading trolley underwheel height	Min 44 ± 1 cm Max 55 ± 1 cm

#### Class I MD compliant with UE Reg. 2017/745

ELLA SELF TF BLUE WHEELS Ø200  
MM

ST10608C

CND Classification V0899

Registration number 159749



Rev.0 (03/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)