

ELLA SELF

Self-loading stretcher



Ella series field stretchers Ella Self are devices for lifting and carrying people for use in field rescue situations. The foldable, lockable legs allow for transferring either by means of wheels on the legs or by lifting when the legs are locked in the closed position. Parking is permitted with the stretcher alone or with a trolley with closed legs.

Specific features	Technical data	
Polyethylene board easy to clean	Length with handles closed	1980 ±10 mm
Color coded commands	Width	570 ± 10 mm
Adjustable backrest	Wheels	Ø 200 mm
Board with Trendelenburg antishock position	Dimensions sidebars	680 x 200 mm
2 fixed wheels and 2 pivoting with large diameter to make movement easier	Materials	Steel/Nylon/Al/PE
Brakes on rear wheels	Weight	40 ± 1 kg
Foldable side bars Safe bar	Capacity	170 kg
System to lock the legs in closed position to lift and carry the device without using the wheels	Loading trolley underwheel height	Min 44 ± 1 cm Max 55 ± 1 cm
Equipped with telescopic handles to make lifting operations easier		

Accessories	
EN90003C	END-T Table for stretchers 10 g certified
ST00491A	STX 90 Telescopic headrest for stretchers
IF01044B	Track 4 telescopic IV holder for Ella Self
CB09025C	Tanker oxygen tank holder for stretcher side bar
CB09028C	Tanker support for Safe Bar
ST70000A	QMX 777 Black PVC anatomic mattress
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Thorax Belt w/integrated re-winding system
ST42100A	FP Posterior fastener for self-loading stretchers
ST42200A	FA Anterior fastener for self-loading stretchers

Standard equipment
STX 592 – Two piece restraints with metal buckle

Class I MD compliant with UE Reg. 2017/745	
ELLA SELF T SILVER TELESCOPIC HANDLES WHEELS °200M	
ST10701C	
CND Classification	V0899
Registration number	159730

