

ELLA SELF

Self-loading stretcher



Ella series field stretchers Ella Self are devices for lifting and carrying people for use in field rescue situations. The foldable, lockable legs allow for transferring either by means of wheels on the legs or by lifting when the legs are locked in the closed position. Parking is permitted with the stretcher alone or with a trolley with closed legs.

Specific features

- Polyethylene board easy to clean
- Color coded commands
- Adjustable backrest in 5 positions in addition to the horizontal
- Board with Trendelenburg antishock and Fowler position
- 2 fixed wheels and 2 pivoting with large diameter to make movement easier
- Brakes on rear wheels
- Foldable side bars *Safe bar*
- System to lock the legs in closed position to lift and carry the device without using the wheels

Accessories

EN90003C	END-T Table for stretchers 10 g certified
ST00491A	STX 90 Telescopic headrest for stretchers
IF01044B	Track 4 telescopic IV holder for Ella Self
CB09025C	Tanker oxygen tank holder for stretcher side bar
CB09028C	Tanker support for Safe Bar
ST70000A	QMX 777 Black PVC anatomic mattress
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Thorax Belt w/integrated re-winding system
ST42100A	FP Posterior fastener for self-loading stretchers
ST42200A	FA Anterior fastener for self-loading stretchers

Technical data

Length with handles closed	1980 ±10 mm
Width	570 ± 10 mm
Wheels	Ø 200 mm
Dimensions sidebars	680 x 200 mm
Materials	Steel/Nylon/Al/PE
Weight	40 ± 1 kg
Capacity	170 kg
Loading trolley underwheel height	Min 44 ± 1 cm Max 55 ± 1 cm

Class I MD compliant with UE Reg. 2017/745

MOD.B005 ELLA TF BLUE
W/ROLLER MAT.SIST.DIS.LEGS
ST10908C



CND Classification

V0899

Standard equipment

STX 592 – Two piece restraints with metal buckle

Rev.0 (03/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it