

Spencer 427

Transport chair with four wheels and armrests



Spencer 427 is a transport chair for the evacuation and transport of patient from the rescue place to the ambulance or to another device intended for transport on the vehicle.

Specific features

- Easy to clean PVC sheets
- 4 rear foldable handles to reduce encumbrance and to facilitate the lifting in every situation
- Telescopic front handles to facilitate the handling
- Foldable to save space
- Strong and light frame
- Foldable footrest to increase stability and comfort of the patient
- Front wheels with brake for rotation and axis

Standard equipment

- | | |
|----------|-----------------------------------|
| ST00423A | Couple of belts for chair 420/425 |
|----------|-----------------------------------|

Accessories

- | | |
|----------|---------------------------------------------------|
| ST30451B | MOUNTING BRACKETS FOR CHAIR SPENCER 450-55/420-25 |
| ST36452E | UPPER FIXATION SYSTEM FOR TRANSPORT CHAIRS |

Technical data

Height	1200 ± 10mm
Width	490 ± 10mm
Length opened w/closed handles	680 ± 10mm
Length opened w/opened handles	1180 ± 10mm
Thickness folded	210 ± 10mm
Rear wheels	ø125 mm
Front wheels	ø100 mm
Materials	Steel/Al/PU/PVC
Weight	11 kg ± 0,3
Loading capacity	150 kg

Class I Med.Device according to Dir. 93/42/CEE

SPENCER 427 SILVER 4 WHEELS
W/ARMRESTS
ST11425A



CND Classification	V9099
Registration number	93519

Rev.3 (02/04/2020)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
[.spencer.it](http://www.spencer.it)