

Spencer 420

Transport chair with two wheels



Spencer 420 is a transport chair for the evacuation and transport of patient from the rescue place to the ambulance or to another device intended for transport on the vehicle.

Specific Characteristics

Easy to clean PVC sheets
 2 rear foldable handles to reduce encumbrance and to facilitate the lifting in every situation
 Telescopic front handles to facilitate the handling
 Foldable to save space
 Strong and light frame

Standard equipment

ST00423A Pair of belts for chair 420/425 orange

Accessories

ST36451A F400 20G FASTENER

Technical data

Length opened w/closed handles	660 ± 10 mm
Length opened w/opened handles	1150 ± 10 mm
Width	510 ± 10 mm
Height opened	900 ± 10 mm
Thickness folded	210 ± 10 mm
Rear wheels	ø125 mm
Materials	Steel/Al/PU/PVC
Weight	10,5 kg ± 0,3
Loading capacity	150 kg

Class I Med.Device according to Regulation 2017/745/UE

STAIR CHAIR SILVER SPENCER 420S

ST20420A

CND Classification



V9099

Rev.0 (02/08/2022)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it